



# MD300K1 Veterinary Pulse Oximeter

- \* Real-time accurately monitoring the data of SPO2 and PR;
- \* Motion interference resistance, stable performance;
- \* All digital blood-oxygen technology , SpO2 and Pulse rate can be shown correctly in low perfusion condition;
- \* Up to 99 IDs codes, can storage 72 hours data and upload to computer to retrieve data using analytic software;
- \* High brightness LED display;
- \* Low power consumption: using 2 AA batteries, nominal duration of service not less than 24 hours;
- \* Visual & audio Alarm: alarming when measured value exceed default range of SPO2 and PR;
- \* Humanized Design: rubber belt make it possible to fix on bedstead.



## Beijing Choice Electronic Technology Co., Ltd

Address: Rm 1127-1128, Bailangyuan Building B, Fuxing Road A36, Beijing, P. R. China Zip: 100039

Tel: +86 10 88204188, 68221006, 88203520 FAX: +86 10 88204632

Email: market@choicemmed.com

Website: www.choicemmed.com